



Care Instructions for Silicone/Gel Liners

A silicone/gel liner should last approximately 9 months.

How to look after your silicone/gel liner:

- 1) Wash your liner daily using mild unscented soap and water solution, gel/silicone side only. Thoroughly rinse; then dry with a lint free cloth/towel.
- 2) Twice a week clean silicone/gel side of liner with an alcohol wipe or diluted antiseptic solution.
- 3) Leave the liner fabric side out to dry fully. Leaving silicone/gel side out can cause liner material to crack and fail.
- 4) Leave liner to dry naturally as any form of direct heat will cause the liner material to fail.

If your liner does not last approximately 9 months, please tell your Prosthetist so any potential issues causing liner failure can be addressed.